

Medical Gas Pipeline

Medical Gas Control Panel

The SILBERMANN medical gas control panel supplies adjustable pressure, as required, to pneumatic-driven surgical tools in the operating theaters.

Silbermann's gas control panel is designed to meet the latest and most rigid American NFPA99 and European EN ISO 7396-1 requirements.





Features

- Compatible with Nitrogen or Industrial Air
- Easy, safe and reliable connection between the medical equipment and the medical gas distribution system
- Anodized Aluminum front with rounded edges for aesthetic appeal
- Manual shut-off valve
- Inlet & outlet pressure gauges
- DISS or NIST outlet connection
- A Knob to adjust outlet pressure

Medical Gas System Standards:

NFPA99, EN ISO 7396-1

Silbermann's Quality Management System:

- ISO 9001-2015
- ISO 13485-2016
- CE mark








Technical Data

- Panel: 3 mm anodized Aluminum plate
- Compatible with Nitrogen or Instrument Air
- Box: Galvanized steel case, powder coated Ral 9016 (White)
- Inlet Pressure Guage: Up to 25 bars
- Adjusted Pressure: 0.5-16 bar.
- Flow Capacity: 70 m³/h @ 10 bar
 90 m³/h @ 16 bar



Dimensions [mm]





Catalogue Numbers*

Type of Gas	Exposed/Recessed Panel	Catalog No.
Air	Recessed	7167400A
Nitrogen	Recessed	7167300

^{*} Catalogue No. are for Control Panels with DISS adapters. NIST adapters are available by customers request.

Distributed By:		
ı		

SILBERMANN Technologies Ltd.

Address: #5 Harakevet st. POB 4605, Petah-Tikva, 4900715, ISRAEL

Tel: +972-3-9309186 | Fax: +972-3-9315721 E-Mail: info@silbermnn.com, www.silbermnn.com

